# LogTag®

Temperature Recorder

## vaxtag®+

### ENGLISH

#### User Instructions

#### Activating the vaxtag®+ recorder

The recorder is delivered in 'sleep mode' with the LCD screen blank. Every **vaxtag®+** must be activated before its first use. Press and hold the ◀ and ▶ buttons together. Continue to hold while **READY** flashes on the display; release the buttons when the flashing stops (after approx. 6 seconds). The recorder will now show **CLOCK ADJ** and is ready for the clock to be set.

#### Setting the Display Clock

Once activated, the display shows **CLOCK ADJ**.

Initially the minute digits flash. Press ◀ to increment and ▶ to store.

Now the hour digits flash. Press ◀ to increment and ▶ to store.

Now the day of the month flashes. Press ◀ to toggle and ▶ to store (month and year are shown for reference).

The clock can be adjusted at any time by pressing and holding the ◀ and ▶ buttons simultaneously. Continue to hold while **CLOCK ADJ** flashes on the display; release the buttons when the flashing stops and adjust the clock. The **vaxtag®+** recorder displays a 24 hour clock.

#### Starting the Recorder

The display must show **READY** before the recorder can be started.

Press and hold the Start ▶ button. Continue to hold while **STARTING** flashes on the display; release the button when the flashing stops (after approx 5 seconds).

The **vaxtag®+** recorder is now counting down a start delay of 30 minutes and is ready to be placed in the area to monitor. At the end of this period the **vaxtag®+** recorder will start collecting data and show **[REC•]** on the display.

Continue to read these instructions to learn about the display and operation of the **vaxtag®+**

#### “Paused” Function

Pressing either button while recording excludes the next five readings from the alarm and statistics calculations; **PAUSED** will now be shown on the display. This allows the user to review the statistics or clear an alarm without causing an invalid reading, alarm or statistic.

#### Recording Display

While recording the display is updated every 6 minutes, showing the most recent temperature reading taken, the current time, battery status and combined alarm status and history.

#### Examples:

Display shows alarms occurring 7 and 19 days ago. Current temperature is over the upper alarm limit of +8°C (indicated by the up-arrow) but the duration is not enough to trigger an alarm (✓).

10 hours later, the over-temperature duration triggers an alarm. The **ALARM** symbol, the cross symbol (X) and the **T** marker are shown.

A further 2 hours later the temperature has returned to the acceptable range but the alarm remains. The day summary has shifted by 1 day as the time has gone through midnight.

#### Clearing an Alarm

When an alarm event occurs, **ALARM** symbol and X are displayed. To clear the **ALARM** symbol and X, press and hold the Clear ▶ button. Alarm events are stored and can be reviewed with the day statistics.

#### Reviewing Day Statistics

Press the Review ◀ button. Each press of the button steps through the daily Max and Min temperatures starting with today, moving backwards through the days (the displayed day will flash).

**T** (flashing) Max Temp (Normal range)

**T** (flashing) Min Temp (Normal range)

Press ◀

**T** (flashing) Max Temp, is above +8°C and has triggered an alarm. Duration was 10:30

**T** (flashing) Min Temp (Normal range)

Press ◀

**T** (flashing) Min Temp (Normal range)

If less than 30 days of recordings have been collected, the day marker only advances to the earliest day of recording. After the statistics for the earliest day have been displayed, the next press of the ◀ button rolls the summary back to the statistic for today **T**.

Pressing the ▶ button at any time shows the normal recording display again.

#### PDF

You can view a PDF of the recorded data by plugging the recorder into a USB socket of any PC capable of displaying PDF files. PDF reader software is required, such as Acrobat Reader or similar.

The PDF contains a summary of the trip, alarm details, a chart, and a list of the recorded temperatures.

It is possible to plug the **vaxtag®+** recorder into a USB socket while it is still recording, but not recommended. You will be able to view the PDF file, but during this time the recorder will not record any temperature data, and **USB** will be shown.

#### Battery

The **vaxtag®+** recorder contains a Lithium Battery. The end of the battery life is indicated by a low battery symbol, the recorder should be replaced within 30 days when this symbol appears.

Dispose or recycle the battery in accordance to your local regulations.

Do not expose the recorder to extreme temperatures as it may lead to the destruction of the battery and may cause injuries. Keep out of reach of children.

#### Liability

The manufacturer shall not be held liable:

- If the device was used beyond the manufacturer's given limitations.
- For any claims due to the improper storage and use of the device.
- For any problems with the cooling unit.
- For the bad quality of the monitored goods, if any.
- For incorrect readings if the device was used with activated low battery sign.

#### Useful life

The operational life of the **vaxtag®+** is approx. 2 years after activation subject to:

- The device was not stored for more than 2 years prior to activation.
- The read out mode is not excessively used e.g. several times a day.
- Storing & operating of the device remains inside the recommendations of the manufacturer.

Attention: The **vaxtag®+** monitors temperature exposure and not the product quality. Its purpose is to signal if product quality evaluation/testing is required.

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